

Compound T129457

Chemical Properties


CAS No. :

Formula: C48H76O18

Molecular Weight: 941.12

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(2S,3R,4S,5R,6R)-6-(((2R,3R,4S,5S,6S)-2-(((3S,6aR,6bS,8aR,9R,12aS,14bR)-9-hydroxy-4,4,6a,6b,8a,11,11,14b-octamethyl-1,2,3,4,4a,5,6,6a,6b,7,8,8a,9,10,11,12,12a,14,14a,14b-icosahydricen-3-yl]oxy)-6-carboxy-4,5-dihydroxyoxan-3-yl]oxy)-3,4-dihydroxy-5-(((2S,3R,4R,5R,6S)-3,4,5-trihydroxy-6-methyloxan-2-yl]oxy)oxane-2-carboxylic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0626 mL	5.3128 mL	10.6256 mL
5 mM	0.2125 mL	1.0626 mL	2.1251 mL
10 mM	0.1063 mL	0.5313 mL	1.0626 mL
50 mM	0.0213 mL	0.1063 mL	0.2125 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481