

Compound T127758

Chemical Properties


CAS No. :

Formula: C₂₈H₂₆O₁₃

Molecular Weight: 570.5

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

(2S,3R,4S,5R,6R)-3,5-dihydroxy-2-[(5-hydroxy-4-oxo-2-phenyl-3,4-dihydro-2H-1-benzopyran-7-yl)oxy]-6-(hydroxymethyl)oxan-4-yl 3,4,5-trihydroxybenzoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7528 mL	8.7642 mL	17.5285 mL
5 mM	0.3506 mL	1.7528 mL	3.5057 mL
10 mM	0.1753 mL	0.8764 mL	1.7528 mL
50 mM	0.0351 mL	0.1753 mL	0.3506 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481