

## Fmoc-Gln(Trt)-Ser(ψ(Me,Me)pro)-OH

## Chemical Properties

CAS No. :	1821378-64-6
Formula:	C <sub>45</sub> H <sub>43</sub> N <sub>3</sub> O <sub>7</sub>
Molecular Weight:	737.839
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Fmoc-Gln(Trt)-Ser(ψ(Me,Me)pro)-OH is a dipeptide.
Targets(IC50)	Amino Acids and Derivatives

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3553 mL	6.7765 mL	13.5531 mL
5 mM	0.2711 mL	1.3553 mL	2.7106 mL
10 mM	0.1355 mL	0.6777 mL	1.3553 mL
50 mM	0.0271 mL	0.1355 mL	0.2711 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481