

## (2E,4E,6E)-octa-2,4,6-trienoic acid

## Chemical Properties

CAS No. :

Formula: C<sub>8</sub>H<sub>10</sub>O<sub>2</sub>

Molecular Weight: 138.17

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	(2E,4E,6E)-octa-2,4,6-trienoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.2375 mL	36.1873 mL	72.3746 mL
5 mM	1.4475 mL	7.2375 mL	14.4749 mL
10 mM	0.7237 mL	3.6187 mL	7.2375 mL
50 mM	0.1447 mL	0.7237 mL	1.4475 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481