

Pomegralignan

Chemical Properties


CAS No. : 1562492-32-3

Formula: C₂₄H₂₈O₁₂

Molecular Weight: 508.47

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Pomegralignan is a natural product for research related to life sciences. The catalog number is TN5430 and the CAS number is 1562492-32-3.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9667 mL	9.8334 mL	19.6668 mL
5 mM	0.3933 mL	1.9667 mL	3.9334 mL
10 mM	0.1967 mL	0.9833 mL	1.9667 mL
50 mM	0.0393 mL	0.1967 mL	0.3933 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481