


Notoginsenoside S

Chemical Properties

CAS No. :	575446-95-6
Formula:	C ₆₃ H ₁₀₆ O ₃₀
Molecular Weight:	1343.51
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Notoginsenoside S may have hepatoprotective effects.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7443 mL	3.7216 mL	7.4432 mL
5 mM	0.1489 mL	0.7443 mL	1.4886 mL
10 mM	0.0744 mL	0.3722 mL	0.7443 mL
50 mM	0.0149 mL	0.0744 mL	0.1489 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481