


Phenyl salicylate

Chemical Properties

CAS No. :	118-55-8
Formula:	C13H10O3
Molecular Weight:	214.22
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Phenyl salicylate (NSC-33406) is a chemical substance. It is synthesized by heating salicylic acid with phenol.
Targets(IC50)	Others,Antibacterial

Solubility Information

Solubility	DMSO: 25 mg/mL (116.7 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.6681 mL	23.3405 mL	46.681 mL
5 mM	0.9336 mL	4.6681 mL	9.3362 mL
10 mM	0.4668 mL	2.334 mL	4.6681 mL
50 mM	0.0934 mL	0.4668 mL	0.9336 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Lapczynski A , Jones L , Mcginty D , et al. Fragrance material review on phenyl salicylate[J]. Food & Chemical Toxicology, 2007, 45(1-suppl-S1):S472-S476.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481