

Compound N021-0003

Chemical Properties


CAS No. :

Formula: C₂₄H₃₂N₂O₇

Molecular Weight: 460.527

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N021-0003 is a useful organic compound for research related to life sciences and the catalog number is T131319.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1714 mL	10.8571 mL	21.7141 mL
5 mM	0.4343 mL	2.1714 mL	4.3428 mL
10 mM	0.2171 mL	1.0857 mL	2.1714 mL
50 mM	0.0434 mL	0.2171 mL	0.4343 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481