


Fmoc-Phe-Lys(Trt)-PAB-PNP

Chemical Properties

CAS No. :	1116086-09-9
Formula:	C63H57N5O9
Molecular Weight:	1028.175
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Fmoc-Phe-Lys(Trt)-PAB-PNP is a cleavable linker designed specifically for antibody-drug conjugate (ADC) synthesis.
Targets(IC50)	ADC Linker,PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9726 mL	4.863 mL	9.7259 mL
5 mM	0.1945 mL	0.9726 mL	1.9452 mL
10 mM	0.0973 mL	0.4863 mL	0.9726 mL
50 mM	0.0195 mL	0.0973 mL	0.1945 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481