


Mepixanox

Chemical Properties

CAS No. :	17854-59-0
Formula:	C ₂₀ H ₂₁ NO ₃
Molecular Weight:	323.39
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Mepixanox is an analeptic drug. It is used for treating respiratory and cardiorespiratory insufficiency.
Targets(IC50)	Others,CRFR

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0922 mL	15.4612 mL	30.9224 mL
5 mM	0.6184 mL	3.0922 mL	6.1845 mL
10 mM	0.3092 mL	1.5461 mL	3.0922 mL
50 mM	0.0618 mL	0.3092 mL	0.6184 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Grossi G, et al. High-performance liquid chromatographic determination of mepixanthone in serum. J Chromatogr. 1984 Jul 13;309(1):214-8.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481