

2-[(2E)-2-(hydroxyimino)-4-methylpentanamido]-4-methylpentanoic acid

Chemical Properties

CAS No. :

Formula: C₁₂H₂₂N₂O₄

Molecular Weight: 258.32

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	2-[(2E)-2-(hydroxyimino)-4-methylpentanamido]-4-methylpentanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8712 mL	19.3558 mL	38.7117 mL
5 mM	0.7742 mL	3.8712 mL	7.7423 mL
10 mM	0.3871 mL	1.9356 mL	3.8712 mL
50 mM	0.0774 mL	0.3871 mL	0.7742 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481