

L 698680

## Chemical Properties


CAS No. : 139548-02-0

Formula: C<sub>17</sub>H<sub>18</sub>N<sub>2</sub>O<sub>3</sub>

Molecular Weight: 298.34

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	L 698680 is a 6-hydroxy metabolite of L-696229.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3519 mL	16.7594 mL	33.5188 mL
5 mM	0.6704 mL	3.3519 mL	6.7038 mL
10 mM	0.3352 mL	1.6759 mL	3.3519 mL
50 mM	0.067 mL	0.3352 mL	0.6704 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Lee LL, Herold ML, Zacchei AG. High-performance liquid chromatographic method for the determination of an HIV-1 non-nucleoside reverse transcriptase inhibitor (L-696,229) in plasma samples from animals. J Chromatogr B Biomed Appl. 1996 Oct 25;685(2):323-8. PubMed PMID: 8953174.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481