

Compound T129068

Chemical Properties

CAS No. :

Formula: C34H54O23

Molecular Weight: 830.78

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	[(2R,3S,4S,5R,6S)-6-(((2R,3R,4S,5S,6R)-2-(((2R,3R,4S,5S,6R)-4,5-dihydroxy-6-(hydroxymethyl)-3-(((2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy)oxan-2-yl]oxy)-4,5-dihydroxy-6-(((2E)-2-methylbut-2-enoyl]oxy)methyl)oxan-3-yl]oxy)-3,4,5-trihydroxyoxan-2-yl]methyl (2E)-2-methylbut-2-enoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2037 mL	6.0184 mL	12.0369 mL
5 mM	0.2407 mL	1.2037 mL	2.4074 mL
10 mM	0.1204 mL	0.6018 mL	1.2037 mL
50 mM	0.0241 mL	0.1204 mL	0.2407 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481