

## Neuraminidase-IN-4

## Chemical Properties

CAS No. :

Formula: C<sub>21</sub>H<sub>20</sub>N<sub>2</sub>O<sub>6</sub>S

Molecular Weight: 428.46

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Neuraminidase-IN-4 (Compound 4b) is a potent inhibitor of neuraminidase (EC <sub>50</sub> : 1.59 $\mu$ M) and exhibits significant antiviral activity against A/chicken/Hubei/327/2004 (H5N1-DW).
Targets(IC <sub>50</sub> )	Others, Influenza Virus

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3339 mL	11.6697 mL	23.3394 mL
5 mM	0.4668 mL	2.3339 mL	4.6679 mL
10 mM	0.2334 mL	1.167 mL	2.3339 mL
50 mM	0.0467 mL	0.2334 mL	0.4668 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481