

2-Chloro-N6-iso-pentenyladenosine

Chemical Properties

CAS No. :


Formula:

Molecular Weight:

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-Chloro-N6-iso-pentenyladenosine is a nucleoside derivative, specifically a halo-nucleoside, characterized by modifications at the 2 and 6 positions of the purine nucleoside.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481