


MMG-0358

Chemical Properties

CAS No. : 1378976-02-3
 Formula: C₈H₆ClN₃O
 Molecular Weight: 195.61
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	MMG-0358 is a novel potent IDO1 inhibitor, showing low cytotoxicity and higher selectivity for IDO1 over TDO enzyme.
Targets(IC50)	Others,IDO,Indoleamine 2,3-Dioxygenase (IDO)

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.1122 mL	25.5611 mL	51.1221 mL
5 mM	1.0224 mL	5.1122 mL	10.2244 mL
10 mM	0.5112 mL	2.5561 mL	5.1122 mL
50 mM	0.1022 mL	0.5112 mL	1.0224 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481