

Benzyl glycinate hydrochloride

Chemical Properties

CAS No. : 2462-31-9


Formula: C₉H₁₂ClNO₂

Molecular Weight: 201.65

Storage: Keep away from moisture

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Benzyl glycinate hydrochloride is a glycine derivative.
Targets(IC50)	Amino Acids and Derivatives

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.9591 mL	24.7954 mL	49.5909 mL
5 mM	0.9918 mL	4.9591 mL	9.9182 mL
10 mM	0.4959 mL	2.4795 mL	4.9591 mL
50 mM	0.0992 mL	0.4959 mL	0.9918 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481