

Compound 0407-0029

## Chemical Properties


CAS No. :

Formula: C<sub>30</sub>H<sub>43</sub>NO<sub>4</sub>

Molecular Weight: 481.677

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Compound 0407-0029 is a useful organic compound for research related to life sciences and the catalog number is T131678.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0761 mL	10.3803 mL	20.7607 mL
5 mM	0.4152 mL	2.0761 mL	4.1521 mL
10 mM	0.2076 mL	1.038 mL	2.0761 mL
50 mM	0.0415 mL	0.2076 mL	0.4152 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481