

TPP-1 Acetate(2426685-25-6 Free base)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Store at low temperature,Keep away from moisture

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

TPP-1 Acetate(2426685-25-6 Free base) is a peptide molecule that can be used in life science related research.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481