


2-(2-Thienyl)-1,3-thiazole-4-carboxylic acid

Chemical Properties

CAS No. :	24044-07-3
Formula:	C ₈ H ₅ NO ₂ S ₂
Molecular Weight:	211.261
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-(2-Thienyl)-1,3-thiazole-4-carboxylic acid, with CAS No. 24044-07-3, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 2-(2-Thienyl)-1,3-thiazole-4-carboxylic acid provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.7335 mL	23.6675 mL	47.335 mL
5 mM	0.9467 mL	4.7335 mL	9.467 mL
10 mM	0.4734 mL	2.3668 mL	4.7335 mL
50 mM	0.0947 mL	0.4734 mL	0.9467 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481