

Degalactotigonin

Chemical Properties


CAS No. : 39941-51-0

Formula: C₅₀H₈₂O₂₂

Molecular Weight: 1035.184

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

Degalactotigonin is a natural product.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.966 mL	4.8301 mL	9.6602 mL
5 mM	0.1932 mL	0.966 mL	1.932 mL
10 mM	0.0966 mL	0.483 mL	0.966 mL
50 mM	0.0193 mL	0.0966 mL	0.1932 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481