


Benzoylmalic acid

Chemical Properties

CAS No. :	22138-51-8
Formula:	C ₁₁ H ₁₀ O ₆
Molecular Weight:	238.195
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Benzoylmalic acid is a useful organic compound for research related to life sciences. The catalog number is T125309 and the CAS number is 22138-51-8.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.1982 mL	20.9908 mL	41.9815 mL
5 mM	0.8396 mL	4.1982 mL	8.3963 mL
10 mM	0.4198 mL	2.0991 mL	4.1982 mL
50 mM	0.084 mL	0.4198 mL	0.8396 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481