


## 3'-O-Methylmurraol

## Chemical Properties

CAS No. : 1891097-17-8  
 Formula: C<sub>16</sub>H<sub>18</sub>O<sub>4</sub>  
 Molecular Weight: 274.31  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



## Biological Description

Description	3'-O-Methylmurraol may have potential anti-inflammatory activity, it exhibits inhibition (IC <sub>50</sub> ≤ 7.31 ug/mL) of superoxide anion generation by human neutrophils in response to formyl-L-methionyl-L-leucyl-L-phenylalanine/cytochalasin B (fMLP/CB).
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6455 mL	18.2276 mL	36.4551 mL
5 mM	0.7291 mL	3.6455 mL	7.291 mL
10 mM	0.3646 mL	1.8228 mL	3.6455 mL
50 mM	0.0729 mL	0.3646 mL	0.7291 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481