


**(±)-Isovalerylshikonin****Chemical Properties**

CAS No. :	52387-14-1
Formula:	C <sub>21</sub> H <sub>24</sub> O <sub>6</sub>
Molecular Weight:	372.417
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.

**Biological Description**

Description	(±)-Isovalerylshikonin is a useful organic compound for research related to life sciences. The catalog number is T131720 and the CAS number is 52387-14-1.
Targets(IC50)	Others,Antibacterial,Antibiotic

**Preparing Stock Solutions**

	1mg	5mg	10mg
1 mM	2.6851 mL	13.4257 mL	26.8514 mL
5 mM	0.537 mL	2.6851 mL	5.3703 mL
10 mM	0.2685 mL	1.3426 mL	2.6851 mL
50 mM	0.0537 mL	0.2685 mL	0.537 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481