

QNN33358

Chemical Properties


CAS No. : 77133-35-8

Formula: C₂₁H₄₂O₄

Molecular Weight: 358.56

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	QNN33358 is an alkylacylglycerol acts as a protein kinase C (PKC) inhibitor and diacylglycerol (DAG) antagonist.
Targets(IC50)	Others,mTOR

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7889 mL	13.9447 mL	27.8893 mL
5 mM	0.5578 mL	2.7889 mL	5.5779 mL
10 mM	0.2789 mL	1.3945 mL	2.7889 mL
50 mM	0.0558 mL	0.2789 mL	0.5578 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481