


## 2-Pyrrolidin-1-ylisonicotinonitrile

## Chemical Properties

CAS No. :	127680-87-9
Formula:	C10H11N3
Molecular Weight:	173.2144
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2-Pyrrolidin-1-ylisonicotinonitrile, with CAS No. 127680-87-9, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 2-Pyrrolidin-1-ylisonicotinonitrile provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.7733 mL	28.8667 mL	57.7334 mL
5 mM	1.1547 mL	5.7733 mL	11.5467 mL
10 mM	0.5773 mL	2.8867 mL	5.7733 mL
50 mM	0.1155 mL	0.5773 mL	1.1547 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481