


3-Hydroxyisovaleric acid

Chemical Properties

CAS No. :	625-08-1
Formula:	C ₅ H ₁₀ O ₃
Molecular Weight:	118.13
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	3-Hydroxyisovaleric acid (Beta-Hydroxyisovaleric acid) is an endogenous metabolite excreted in the urine and can serve as an early and sensitive indicator of biotin deficiency.
Targets(IC50)	Endogenous Metabolite

Solubility Information

Solubility	DMSO: 50 mg/mL (423.26 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (16.93 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	8.4653 mL	42.3263 mL	84.6525 mL
5 mM	1.6931 mL	8.4653 mL	16.9305 mL
10 mM	0.8465 mL	4.2326 mL	8.4653 mL
50 mM	0.1693 mL	0.8465 mL	1.6931 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Sabourin PJ, Bieber LL. Arch Biochem Biophys. 1981 Jan;206(1):132-44.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481