


## Azido-PEG10-Boc

## Chemical Properties

CAS No. :

Formula: C27H53N3O12

Molecular Weight: 611.72



Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Azido-PEG10-Boc is a PEG-based linker for PROTACs, essential for combining two key ligands necessary for the formation of PROTAC molecules, enabling selective protein degradation through the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two ligands linked together: one binding to an E3 ubiquitin ligase and the other to the target protein. They leverage the intracellular ubiquitin-proteasome system to selectively degrade target proteins [1].

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6347 mL	8.1737 mL	16.3473 mL
5 mM	0.3269 mL	1.6347 mL	3.2695 mL
10 mM	0.1635 mL	0.8174 mL	1.6347 mL
50 mM	0.0327 mL	0.1635 mL	0.3269 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481