


Citrazinic acid

Chemical Properties

CAS No. :	99-11-6
Formula:	C ₆ H ₅ NO ₄
Molecular Weight:	155.11
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Citrazinic acid could inhibit 3-dehydroquinate dehydratase from Mycobacterium tuberculosis. It was also used to synthesize the derivatives of pyrimidinones and oxazinones.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.447 mL	32.2352 mL	64.4704 mL
5 mM	1.2894 mL	6.447 mL	12.8941 mL
10 mM	0.6447 mL	3.2235 mL	6.447 mL
50 mM	0.1289 mL	0.6447 mL	1.2894 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Sarkar S, Chowdhury J, Dutta S, Pal T. A pH dependent Raman and surface enhanced Raman spectroscopic studies of citrazinic acid aided by theoretical calculations. *Spectrochim Acta A Mol Biomol Spectrosc.* 2016 Dec 5;169:108-15. doi: 10.1016/j.saa.2016.06.023. Epub 2016 Jun 16. PubMed PMID: 27344521.

Hossan AS, Abu-Melha HM, Al-Omar MA, Amr Ael-G. Synthesis and antimicrobial activity of some new pyrimidinone and oxazinone derivatives fused with thiophene rings using 2-chloro-6-ethoxy-4-acetylpyridine as starting material. *Molecules.* 2012 Nov 19;17(11):13642-55. doi: 10.3390/molecules171113642. PubMed PMID: 23165308.

Masoud MS, Ibrahim AA, Khalil EA, El-Marghany A. Spectral properties of some metal complexes derived from uracil-thiouracil and citrazinic acid compounds. *Spectrochim Acta A Mol Biomol Spectrosc.* 2007 Jul;67(3-4):662-8. Epub 2006 Aug 25. PubMed PMID: 17029947.

Dias MV, Snee WC, Bromfield KM, Payne RJ, Palaninathan SK, Ciulli A, Howard NI, Abell C, Sacchettini JC, Blundell TL. Structural investigation of inhibitor designs targeting 3-dehydroquinase dehydratase from the shikimate pathway of *Mycobacterium tuberculosis*. *Biochem J.* 2011 Jun 15;436(3):729-39. doi: 10.1042/BJ20110002. PubMed PMID: 21410435.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481