

DBCO-PEG4-alcohol

Chemical Properties

CAS No. : 1416711-60-8


Formula: C₂₉H₃₆N₂O₆

Molecular Weight: 508.61

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	DBCO-PEG4-alcohol is a PEG-based linker for PROTACs that joins two essential ligands, crucial for forming PROTAC molecules, enabling selective protein degradation by leveraging the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs, consisting of two distinct ligands joined by a linker—one for an E3 ubiquitin ligase and the other for the target protein—utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9661 mL	9.8307 mL	19.6614 mL
5 mM	0.3932 mL	1.9661 mL	3.9323 mL
10 mM	0.1966 mL	0.9831 mL	1.9661 mL
50 mM	0.0393 mL	0.1966 mL	0.3932 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481