


Drimenol

Chemical Properties

CAS No. :	468-68-8
Formula:	C ₁₅ H ₂₆ O
Molecular Weight:	222.37
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Drimenol shows antifungal activity, it is able to affect Botrytis cinerea growth with the EC50 value of 80 ppm. Drimendiol also shows quorum sensing inhibition activity.
Targets(IC50)	Antifection

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.497 mL	22.485 mL	44.9701 mL
5 mM	0.8994 mL	4.497 mL	8.994 mL
10 mM	0.4497 mL	2.2485 mL	4.497 mL
50 mM	0.0899 mL	0.4497 mL	0.8994 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481