


Teucardoside

Chemical Properties

CAS No. :	78280-92-9
Formula:	C ₂₁ H ₃₀ O ₁₃
Molecular Weight:	490.458
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Teucardoside is a natural product that can be used as a reference standard. The CAS number of Teucardoside is 78280-92-9.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0389 mL	10.1945 mL	20.389 mL
5 mM	0.4078 mL	2.0389 mL	4.0778 mL
10 mM	0.2039 mL	1.0195 mL	2.0389 mL
50 mM	0.0408 mL	0.2039 mL	0.4078 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481