

Compound N064-0033

## Chemical Properties


CAS No. :

Formula: C<sub>31</sub>H<sub>50</sub>O<sub>3</sub>

Molecular Weight: 470.738

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Compound N064-0033 is a useful organic compound for research related to life sciences and the catalog number is T131207.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1243 mL	10.6216 mL	21.2431 mL
5 mM	0.4249 mL	2.1243 mL	4.2486 mL
10 mM	0.2124 mL	1.0622 mL	2.1243 mL
50 mM	0.0425 mL	0.2124 mL	0.4249 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481