


A 77-01

Chemical Properties

CAS No. : 607737-87-1
Formula: C₁₈H₁₄N₄
Molecular Weight: 286.33
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	A 77-01 is a potent inhibitor of TGF-(beta) type I receptor superfamily activin-like kinase ALK5 with IC ₅₀ of 25 nM.
Targets(IC ₅₀)	ALK,TGF-beta/Smad

Solubility Information

Solubility	DMSO: 11 mg/mL (38.42 mM),Sonification is recommended <1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4925 mL	17.4624 mL	34.9247 mL
5 mM	0.6985 mL	3.4925 mL	6.9849 mL
10 mM	0.3492 mL	1.7462 mL	3.4925 mL
50 mM	0.0698 mL	0.3492 mL	0.6985 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Tojo M, et al. Y Sci. 2005 Nov;96(11):791-800.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481