


3-Acetylcoumarin

Chemical Properties

CAS No. :	3949-36-8
Formula:	C ₁₁ H ₈ O ₃
Molecular Weight:	188.18
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	3-Acetylcoumarin as a practical ligand for copper-catalyzed CN coupling reactions at room temperature. It has neuroprotective and acaricidal properties.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 25 mg/mL (132.85 mM), (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.3141 mL	26.5703 mL	53.1406 mL
5 mM	1.0628 mL	5.3141 mL	10.6281 mL
10 mM	0.5314 mL	2.657 mL	5.3141 mL
50 mM	0.1063 mL	0.5314 mL	1.0628 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Bingchuan Z , Jinxiang L , Ting L , et al. Study on acaricidal bioactivity and quantitative structure activity relationship of coumarin compounds against *Tetranychus cinnabarinus* Bois.(Acari: Tetranychidae)[J]. Chinese

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481