


Methylsyringol

Chemical Properties

CAS No. :	6638-05-7
Formula:	C ₉ H ₁₂ O ₃
Molecular Weight:	168.19
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Methylsyringol is a natural product from Pinus yunnanensis.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 60 mg/mL (356.74 mM), (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.9457 mL	29.7283 mL	59.4566 mL
5 mM	1.1891 mL	5.9457 mL	11.8913 mL
10 mM	0.5946 mL	2.9728 mL	5.9457 mL
50 mM	0.1189 mL	0.5946 mL	1.1891 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481