


Procainamide hydrochloride

Chemical Properties

CAS No. :	614-39-1
Formula:	C ₁₃ H ₂₁ N ₃ O·HCl
Molecular Weight:	271.79
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Procainamide Hydrochloride (Procapan) is the hydrochloride salt form of procainamide, an amide derivative exhibiting class 1A antiarrhythmic property and analog of procaine. Procainamide hydrochloride reversibly binds to and blocks activated (open) voltage-gated sodium channels, thereby blocks the influx of sodium ions into the cell, which leads to an increase in threshold for excitation and inhibit depolarization during phase 0 of the action potential.
Targets(IC50)	DNA Methyltransferase,Sodium Channel,AChR,Autophagy

Solubility Information

Solubility	DMSO: 45 mg/mL (165.57 mM), H2O: 183.97mM, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6793 mL	18.3966 mL	36.7931 mL
5 mM	0.7359 mL	3.6793 mL	7.3586 mL
10 mM	0.3679 mL	1.8397 mL	3.6793 mL
50 mM	0.0736 mL	0.3679 mL	0.7359 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Wojnarowska Z, et al. J Chem Phys. 2012 Apr 28;136(16):164507.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481