


L-Isoleucine

Chemical Properties

CAS No. :	73-32-5
Formula:	C ₆ H ₁₃ NO ₂
Molecular Weight:	131.17
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	L-Isoleucine (Ile) is an essential branched-chain aliphatic amino acid found in many proteins. It is an isomer of LEUCINE. It is important in hemoglobin synthesis and regulation of blood sugar and energy levels.
Targets(IC50)	Amino Acids and Derivatives,Endogenous Metabolite

Solubility Information

Solubility	H ₂ O: 10 mM, DMSO: Insoluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.6237 mL	38.1185 mL	76.2369 mL
5 mM	1.5247 mL	7.6237 mL	15.2474 mL
10 mM	0.7624 mL	3.8118 mL	7.6237 mL
50 mM	0.1525 mL	0.7624 mL	1.5247 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Korman SH, et al. Clin Chem. 2005 Mar;51(3):610-7.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481