


Dihydroferulic acid

Chemical Properties

CAS No. :	1135-23-5
Formula:	C ₁₀ H ₁₂ O ₄
Molecular Weight:	196.2
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Dihydroferulic acid is used as the intermediates of perfume and essence.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 60 mg/mL (305.81 mM), (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.0968 mL	25.4842 mL	50.9684 mL
5 mM	1.0194 mL	5.0968 mL	10.1937 mL
10 mM	0.5097 mL	2.5484 mL	5.0968 mL
50 mM	0.1019 mL	0.5097 mL	1.0194 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Garcia-Conesa MT, et al. Redox Rep. 1997 Oct-Dec;3(5-6):319-23.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481