


Diphenyl oxide

Chemical Properties

CAS No. : 101-84-8
Formula: C₁₂H₁₀O
Molecular Weight: 170.21
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Diphenyl oxide (Oxydibenzene) is used in the manufacture of spices and dyes.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 10 mM, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.8751 mL	29.3755 mL	58.751 mL
5 mM	1.175 mL	5.8751 mL	11.7502 mL
10 mM	0.5875 mL	2.9375 mL	5.8751 mL
50 mM	0.1175 mL	0.5875 mL	1.175 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Garcia-Domenech R, et al. Mol Divers. 2015 May;19(2):357-66.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481