


1-Stearoyl-2-linoleoyl-sn-glycero-3-phosphocholine

Chemical Properties

CAS No. : 27098-24-4
Formula: C44H84NO8P
Molecular Weight: 786.11
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	1-Stearoyl-2-linoleoyl-sn-glycero-3-phosphocholine is an endogenous metabolite.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2721 mL	6.3604 mL	12.7209 mL
5 mM	0.2544 mL	1.2721 mL	2.5442 mL
10 mM	0.1272 mL	0.636 mL	1.2721 mL
50 mM	0.0254 mL	0.1272 mL	0.2544 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481