


AR-13324 M1 metabolite

Chemical Properties

CAS No. :	2309668-15-1
Formula:	C19H19N3O2
Molecular Weight:	321.37
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	AR-13324 M1 metabolite is a hydrolysis metabolite of AR-13324 mesylate.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 41.83 mg/mL (130.16 mM), (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1117 mL	15.5584 mL	31.1168 mL
5 mM	0.6223 mL	3.1117 mL	6.2234 mL
10 mM	0.3112 mL	1.5558 mL	3.1117 mL
50 mM	0.0622 mL	0.3112 mL	0.6223 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481