


PPAR α -MO-1

Chemical Properties

CAS No. : 810677-36-2
Formula: C₂₇H₃₇NO₅
Molecular Weight: 455.59
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	PPAR α -MO-1 is a potent modulator of PPAR α .
Targets(IC50)	PPAR

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.195 mL	10.9748 mL	21.9496 mL
5 mM	0.439 mL	2.195 mL	4.3899 mL
10 mM	0.2195 mL	1.0975 mL	2.195 mL
50 mM	0.0439 mL	0.2195 mL	0.439 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

WO/2004/110982A1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481