


Brexpiprazole-d8

Chemical Properties

CAS No. :	1427049-21-5
Formula:	C ₂₅ H ₂₇ N ₃ O ₂ S
Molecular Weight:	441.62
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Brexpiprazole D8 (OPC-34712 D8) is a deuterium-labeled Brexpiprazole. Brexpiprazole is a partial agonist of human 5-HT _{1A} and dopamine receptors (K _i : 0.12 nM and 0.3 nM).
Targets(IC ₅₀)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2644 mL	11.322 mL	22.6439 mL
5 mM	0.4529 mL	2.2644 mL	4.5288 mL
10 mM	0.2264 mL	1.1322 mL	2.2644 mL
50 mM	0.0453 mL	0.2264 mL	0.4529 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Ishima T, et al. Potentiation of neurite outgrowth by brexpiprazole, a novel serotonin-dopamine activity modulator: a role for serotonin 5-HT_{1A} and 5-HT_{2A} receptors. Eur Neuropsychopharmacol. 2015 Apr;25(4):505-

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481