


GPR40 Activator 1

Chemical Properties

CAS No. : 1309435-60-6
Formula: C₃₁H₃₁NO₃S
Molecular Weight: 497.65
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	GPR40 Activator 1 is a potent GPR40 activator for treatment of type 2 diabetes.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0094 mL	10.0472 mL	20.0944 mL
5 mM	0.4019 mL	2.0094 mL	4.0189 mL
10 mM	0.2009 mL	1.0047 mL	2.0094 mL
50 mM	0.0402 mL	0.2009 mL	0.4019 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Novel spirocyclic compounds , Patent WO 2011066183 A1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481