


IL-17A antagonist 3

Chemical Properties

CAS No. : 2230780-65-9
Formula: C33H33ClN6O4
Molecular Weight: 613.11
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	IL-17A antagonist 3 is an IL-17A antagonist.
Targets(IC50)	IL Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.631 mL	8.1551 mL	16.3103 mL
5 mM	0.3262 mL	1.631 mL	3.2621 mL
10 mM	0.1631 mL	0.8155 mL	1.631 mL
50 mM	0.0326 mL	0.1631 mL	0.3262 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481