


## Kv3 modulator 2

## Chemical Properties

CAS No. :	2101321-76-8
Formula:	C <sub>21</sub> H <sub>21</sub> N <sub>3</sub> O <sub>4</sub>
Molecular Weight:	379.41
Appearance:	no data available
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year



## Biological Description

Description	Kv3 modulator 2 is a potent Kv3 channels modulator.has analgesic activity and is used in the prophylaxis or treatment of related disorders.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6357 mL	13.1784 mL	26.3567 mL
5 mM	0.5271 mL	2.6357 mL	5.2713 mL
10 mM	0.2636 mL	1.3178 mL	2.6357 mL
50 mM	0.0527 mL	0.2636 mL	0.5271 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

## Reference

Giuseppe Alvaro, et al. Hydantoin modulators of kv3 channels.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481