


Mebeverine acid D5

Chemical Properties

CAS No. : 2070015-30-2
Formula: C₁₆H₂₅NO₃
Molecular Weight: 284.41
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Mebeverine acid D5 is the deuterium labeled Mebeverine Acid.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5161 mL	17.5803 mL	35.1605 mL
5 mM	0.7032 mL	3.5161 mL	7.0321 mL
10 mM	0.3516 mL	1.758 mL	3.5161 mL
50 mM	0.0703 mL	0.3516 mL	0.7032 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481