


Mebeverine alcohol

Chemical Properties

CAS No. :	14367-47-6
Formula:	C ₁₆ H ₂₇ NO ₂
Molecular Weight:	265.39
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Mebeverine alcohol is a metabolite of Mebeverine.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 100 mg/mL (376.80 mM), (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.768 mL	18.8402 mL	37.6804 mL
5 mM	0.7536 mL	3.768 mL	7.5361 mL
10 mM	0.3768 mL	1.884 mL	3.768 mL
50 mM	0.0754 mL	0.3768 mL	0.7536 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481