

Compound N039-0043

Chemical Properties


CAS No. :

Formula: C₂₃H₃₄O₃

Molecular Weight: 358.53

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Geminin, Potency: 290.9 nM
-------------	----------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7892 mL	13.9458 mL	27.8917 mL
5 mM	0.5578 mL	2.7892 mL	5.5783 mL
10 mM	0.2789 mL	1.3946 mL	2.7892 mL
50 mM	0.0558 mL	0.2789 mL	0.5578 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481